



CortexCLMS

Calibration Laboratory Management System

Complete User Manual

Covering all modules · Step-by-step guides · Role-based access · Admin configuration



+91 8329925318

✉ cortexaitechnologies@zohomail.in



Pune, Maharashtra, India 411034

Version 1.0 | 2026 | All Rights Reserved

Table of Contents

GETTING STARTED

01	System Overview & Login	4
02	Dashboard	6

CALIBRATION OPERATIONS

03	Jobs — Calibration Work Orders	8
04	Certificates	12
05	Instruments	14
06	Reference Standards	15
07	Calibration Masters	16
08	Environmental Monitoring	17

CRM & BILLING

09	Customers	18
10	Leads & Sales Pipeline	20
11	CRM Activities	21
12	Quotations	22
13	Invoices	23
14	Purchase Orders	24
15	Delivery Challans	25

QUALITY & COMPLIANCE

16	NCR / CAPA (Quality)	26
17	Complaints & Feedback	27
18	Lab Documents	28

19

Internal Audit

29

ADMINISTRATION

20

Tasks

30

21

Engineers

31

22

Users & Role Permissions

32

23

Inventory

34

24

Reports

35

25

Notifications

36

26

Admin Settings

37

27

Customer Portal

40

System Overview & Login

Understanding CLMS architecture, roles, and first-time access

| What is CLMS?

CLMS (Calibration Laboratory Management System) is an integrated web-based platform designed for calibration laboratories. It covers the complete lifecycle — from receiving an instrument for calibration, assigning engineers, issuing certificates, billing customers, managing quality documents, and providing a self-service customer portal.



Calibration Operations

End-to-end job management, engineer assignment, datasheet capture, certificate generation



CRM & Billing

Leads pipeline, quotations, invoices, purchase orders, delivery challans



Quality & Compliance

NCR/CAPA tracking, complaints, document control, internal audits



Customer Portal

Customers can view jobs, download certificates, and raise complaints online

| Login Screen



CLMS — Calibration Laboratory Management System



Sign In to CLMS

Enter your credentials to continue

Email Address

admin@lab.com

Password

.....

Sign In

Customer? [Use the Customer Portal](#) →

User Roles

Every user is assigned one of the following roles. Roles control which modules are visible and what actions are permitted.

Role	Description	Key Access
LAB_ADMIN	Full system administrator	All modules, settings, user management
TECHNICAL_MANAGER	Technical head	Jobs, certificates, engineers, reports
QUALITY_MANAGER	Quality & compliance head	NCR/CAPA, complaints, audits, documents
CALIBRATION_ENGINEER	Performs calibrations	Assigned jobs, datasheets, certificates
SERVICE_ENGINEER	Field/onsite engineer	Onsite jobs, instruments
DATA_ENTRY_OPERATOR	Data entry staff	Customers, instruments, basic data


Note: The LAB_ADMIN sets up the lab during initial registration and cannot be deleted. Additional users are created from *Administration* → *Users*.

Dashboard

At-a-glance KPIs, alerts, and quick navigation

Dashboard Overview

The Dashboard is the first screen after login. It shows real-time KPIs, upcoming due dates, overdue items, and recent activity so the lab manager can take immediate action without navigating through menus.



Dashboard Widgets

Widget	What it Shows	Who Sees It
Active Jobs	Count of jobs not in CLOSED/DELIVERED state	All roles
Pending Review	Jobs in PENDING_REVIEW status awaiting Technical Manager	Admin, Tech Mgr
Certificates Issued	Certificates generated this calendar month	All roles
Overdue Jobs	Jobs past expected delivery date	Admin, Tech Mgr
Status Bar Chart	Job count grouped by current status	Admin, Tech Mgr
Recent Activity Feed	Last 10 events across all modules	All roles

Upcoming Calibrations	Instruments due for re-calibration in next 30 days	Admin, Tech Mgr
Environmental Alerts	Temperature/humidity readings outside acceptable range	Admin, Quality Mgr

CHAPTER 03

Jobs — Calibration Work Orders

Creating, tracking, and managing calibration work orders from intake to delivery

Core Module

Job Lifecycle

Every calibration request becomes a **Job**. A job tracks an instrument through the complete calibration process. Status transitions are strictly controlled to enforce quality workflow.



If corrections are needed after review, the status moves to **CORRECTION_REQUIRED**, and the engineer fixes the datasheet before resubmitting.

Jobs List Screen

Calibration Jobs

+ New Job

+ Bulk Intake

Search by job number, custom

All Statuses

All Engineers

Job No.	Customer / Instrument	Engineer	Status	Actions
JOB-2026-00041	Acme Labs Digital Multimeter	Ravi K.	IN_CALIBRATION	  
JOB-2026-00040	Metro Pharma Weighing Scale	Sanjay M.	APPROVED	  
JOB-2026-00039 BATCH	Tech Solutions Vernier Caliper Onsite	Unassigned	RECEIVED	  

Creating a New Job

Click + **New Job** in the top-right. Fill the form:

Customer	Select the customer from the dropdown (must be added in Customers first)	Required
Instrument	Select the customer's instrument registered in the system	Required
Calibration Procedure	Choose from the grouped procedure library (discipline → sub-category → procedure)	Required
Certificate Type	NABL (accredited) or Non-NABL (in-house). Determines certificate template used	Required
Range / Unit	Auto-filled from procedure; can be overridden per job	
Is Onsite?	Toggle if calibration is performed at the customer's premises	
Visit Date	Required when Is Onsite is enabled	Conditional
PO Number / Challan No.	Customer's purchase order / inward challan reference	
Condition of Item	Describe the item's physical condition at intake	

Bulk Job Intake

Click **+ Bulk Intake** to create multiple jobs under a single batch number for one customer. Select the customer, then multi-select instruments. A single batch covers all jobs created together, making it easy to track related instruments.

Job Actions

Action Button	What it Does	Who Can Use
 Open Job	Opens the Job Workspace (full detail, datasheets, certificate)	All
 Assign Engineer	Opens dialog to assign/re-assign an engineer	Admin, Tech Mgr
 Update Status	Moves the job to the next valid status	Role-dependent
 Delete Job	Permanently deletes the job (with confirmation)	Admin, Tech Mgr

Job Workspace

Click  on any job to open its full workspace. This screen has multiple sections:

- **Job Details** — All metadata: customer, instrument, engineer, dates, remarks
- **Status Progress** — Visual step-by-step pipeline showing current position
- **Datasheet** — Engineers fill in actual calibration readings here
- **Certificate** — Once approved, generate, preview, and download the certificate
- **Audit Trail** — Every status change with timestamp and user

Tip: Engineers can only update status for their own assigned jobs. The Technical Manager or Admin must approve jobs in PENDING_REVIEW status.

CHAPTER 04

Certificates

Viewing, downloading, and verifying calibration certificates

Certificates Module

The Certificates module provides a central registry of all issued calibration certificates. Certificates are automatically generated from approved jobs.

Certificates						
 Search certificate number... All Types ▾ Date Range ▾						
Cert No.	Job No.	Customer	Instrument	Type	Issued Date	Actions
CERT-2026-0041	JOB-2026-00041	Acme Labs	Digital Multimeter	NABL	25 Jun 2026	 

Certificate Types

Type	Template Used	NABL Logo	Use Case
NABL	NABL-accredited template	 Printed	When lab is NABL-accredited for that parameter
Non-NABL	In-house template	 Not printed	Parameters outside NABL scope or interim certs

Certificate Actions

-  **Download** — Downloads the certificate as PDF
-  **Preview** — Opens a print-ready view in a new browser tab
- QR Verification** — Each certificate has a unique QR code. Scanning opens the public verification page confirming authenticity

Public Certificate Verification

Anyone can verify a certificate's authenticity by visiting `/verify/:certNumber` — no login required. This page shows the certificate number, issue date, lab name, and instrument details.

Certificate Number Format: Automatically generated as CERT - {YEAR} - {SEQUENCE}. The format can be customized in Admin Settings → Certificate Settings.

Instruments

Managing customer instruments registered for calibration

| Instruments Registry

The Instruments module stores all customer equipment submitted to the lab for calibration. Each instrument belongs to a customer and has its own calibration history.

| Instrument Fields

Customer	The owner of this instrument	Required
Name	Instrument name (e.g., Digital Multimeter, Vernier Caliper)	Required
Make / Model	Manufacturer and model number	
Serial Number	Instrument's serial number for unique identification	
Range & Unit	Measurement range and unit (e.g., 0–100°C)	
Due Date	Next calibration due date — triggers reminders on dashboard	

Bulk Import: Click *Import CSV* to upload multiple instruments at once using the downloadable template.

| Calibration History

Click the History icon on any instrument to see all past jobs and certificates associated with that instrument.

Reference Standards

Lab's own reference / master equipment calibration tracking

Reference Standards

Reference Standards are the lab's own calibrated equipment used as reference instruments during calibration. They must themselves be calibrated at regular intervals to maintain traceability.

Fields

Name	Name of the reference standard (e.g., Precision Thermometer)	Required
ID / Serial	Lab's internal ID and manufacturer serial	Required
Certificate No.	Calibration certificate number from external calibration	
Valid Until	Expiry date of its current calibration certificate	Required
Traceable To	National/international standard the calibration traces back to	

Alert: Reference standards expiring within 30 days appear highlighted on the Dashboard and in the Notifications panel.

Calibration Masters

Standard procedure library and calibration parameter definitions

Calibration Masters

The Calibration Masters module is the library of calibration procedures the lab follows. These define measurands, units, ranges, and uncertainty parameters used in datasheets and certificates.

Master Fields

Discipline	e.g., Electrical, Mechanical, Thermal, Dimensional	Required
Measurand	The physical quantity being measured (e.g., DC Voltage)	Required
Range	Measurement range (e.g., 0–200 V)	Required
Unit	Unit of measurement (V, mA, °C, mm, etc.)	Required
Uncertainty	Calibration uncertainty value for this range	
Procedure Ref.	Internal/external procedure document reference number	

Environmental Monitoring

Recording temperature, humidity, and lab condition logs

Environmental Monitoring

NABL-accredited labs must maintain environmental records for their calibration rooms. This module allows daily logging of temperature and humidity readings with timestamps.

Log Entry Fields

Date & Time	Timestamp of the reading (auto-filled with current time)	Required
Room / Location	Lab room or area where reading was taken	Required
Temperature (°C)	Dry-bulb temperature reading	Required
Humidity (%RH)	Relative humidity reading	Required
Recorded By	Staff member who took the reading	Required
Remarks	Any abnormalities observed	

Readings outside the acceptable range (typically $23 \pm 2^{\circ}\text{C}$ and 45–65% RH) are flagged in red and generate a dashboard alert.

CHAPTER 09

Customers

Customer database, portal access, and activity timeline

Customer Management

The Customers module is the master directory of all clients. Customer records link to their instruments, jobs, quotations, invoices, and portal access.

Customers						Import CSV	+ New Customer
Code	Name	GSTIN	Email	Phone	Actions		
CUST-001	Acme Laboratories	27AAAAA0000A1Z5	contact@acme.com	+91 9876543210	   		

Customer Fields

Customer Code	Unique code (e.g., CUST-001). Auto-generated or manually set. Cannot be changed after creation.	Required
Name	Company or individual name	Required
GSTIN	GST identification number (used on invoices)	
Email	Primary email — also used for Customer Portal login	
Phone	Contact number	
Address	City, State	

Customer Actions

Button	Action
 Timeline	Opens a drawer showing the customer's full history: quotations received, jobs placed, invoices, payments, and complaints in chronological order
 Set Portal Password	Opens a dialog to set or reset the password the customer uses to log in to the Customer Portal. The customer uses their email + this password.
 Edit	Edit all customer fields (except Code)

 Delete

Permanently delete the customer (blocked if active jobs/invoices exist)

Bulk Import: Use *Import CSV* with the downloadable template to add many customers at once.
Required columns: Customer Code, Name.

CHAPTER 10

Leads & Sales Pipeline

Managing prospective customers from first contact to conversion

| Leads Module

Track potential customers through the sales pipeline before they become formal customers. Each lead moves through stages until won or lost.



| Lead Fields

Company Name	Prospective company name	Required
Contact Person	Name of the decision maker	
Email / Phone	Contact details	
Stage	Current pipeline stage	Required
Source	Where the lead came from: Referral, Cold Call, Website, Exhibition, Other	
Value (₹)	Estimated deal value	
Expected Close	Target date to close the deal	
Notes	Free-text notes about the lead	

CHAPTER 11

CRM Activities

Logging customer interactions and follow-ups

| CRM Activities

Log every customer-facing activity — calls, emails, meetings, and tasks — linked to a specific customer or lead. This creates a complete interaction history.

| Activity Types



| Fields

Type	Activity type from the list above	Required
Customer / Lead	Link to a customer or lead record	Required
Date & Time	When the activity occurred	Required
Subject	Brief title of the interaction	Required
Notes	Detailed notes from the interaction	
Follow-up Date	Next action date	

CHAPTER 12

Quotations

Creating and managing calibration service quotations

| Quotations Module

Generate professional quotations for customers. A quotation can list multiple line items (instrument types and quantities) with rates and GST.

| Quotation Fields

Quotation No.	Auto-generated sequence number	Auto
Customer	Select the customer to quote	Required
Valid Until	Expiry date of this quotation	
Line Items	Description, quantity, unit rate, GST% — multiple rows allowed	Required
Terms & Conditions	Payment terms, delivery scope	

| Quotation Actions

- **Preview / Print** — Formatted A4 PDF with lab letterhead, logo, and signatory
- **Convert to Invoice** — Once accepted, convert directly to an invoice
- **Edit / Delete** — Full CRUD available

CHAPTER 13

Invoices

Billing, GST invoices, and payment tracking

| Invoices Module

Generate GST-compliant tax invoices for calibration services. The system tracks payment status and outstanding dues.

| Invoice Fields

Invoice No.	Auto-generated (e.g., INV-2026-00041)	Auto
Customer	Bill-to customer	Required
Invoice Date	Date of issue	Required
Due Date	Payment due date	
Line Items	Service description, quantity, rate, GST rate (0/5/12/18%)	Required
CGST / SGST / IGST	Auto-calculated based on GST rate and whether customer is intra/inter-state	Auto
Payment Status	UNPAID → PARTIAL → PAID	Auto

Bank account details (account number, IFSC, branch) configured in Admin Settings automatically appear on the invoice PDF.

Purchase Orders

Managing lab procurement and supplier POs

Purchase Orders

Track all purchase orders raised by the lab for consumables, equipment, and services from suppliers.

PO Fields

PO Number	Auto-generated PO reference number	Auto
Supplier	Select the supplier from the customer/vendor directory	Required
PO Date	Date the order was raised	Required
Expected Delivery	Target delivery date	
Line Items	Item description, quantity, unit price	Required
Status	DRAFT → SENT → PARTIAL_RECEIVED → RECEIVED → CANCELLED	

Delivery Challans

Outward and inward goods movement documentation

Delivery Challans

A Delivery Challan (DC) documents the movement of instruments into and out of the lab. It serves as proof of receipt and delivery.

Challan Types

Type	Direction	Use Case
INWARD	Customer → Lab	Customer sends instruments to lab for calibration
OUTWARD	Lab → Customer	Lab returns calibrated instruments to customer

Challan Fields

DC Number	Auto-generated reference number	Auto
Type	INWARD or OUTWARD	Required
Customer	Customer from/to whom instruments are moving	Required
Date	Date of movement	Required
Items	List of instruments with serial numbers and condition	Required
Remarks	Any special instructions or observations	

DC numbers are linked to jobs — the job form has a *Challan No.* field to reference the inward challan when creating a calibration job.

NCR / CAPA — Quality Module

Non-Conformance Reports and Corrective & Preventive Actions

ISO 17025 Compliance

NCR / CAPA

The Quality module manages Non-Conformance Reports (NCR) and Corrective/Preventive Actions (CAPA) as required by ISO/IEC 17025. Any deviation from the defined quality system is documented here.



NCR Fields

NCR Number	Auto-generated (e.g., NCR-2026-0001)	Auto
Title	Brief description of the non-conformance	Required
Category	e.g., Equipment, Method, Personnel, Environment	Required
Description	Detailed description of what went wrong	Required
Root Cause	Root cause analysis (filled during investigation)	
Corrective Action	Steps taken to correct the issue	
Preventive Action	Steps to prevent recurrence	
Target Date	Deadline for CAPA implementation	
Responsible Person	Staff member accountable for closure	

Complaints & Customer Feedback

Managing customer complaints and satisfaction feedback

Complaints Module

All customer complaints — whether raised internally by lab staff or submitted by customers through the Customer Portal — appear here. Each complaint is tracked to closure.

OPEN**IN_PROGRESS****RESOLVED****CLOSED**

Complaints submitted by customers from the Customer Portal automatically appear in this module with the prefix PORTAL/ in the Complaint Number and [Customer: Name] in the description.

Complaint Fields

Complaint No.	Auto-generated (lab: CMP-2026-001, portal: PORTAL/2026/XXXXXXXX)	Auto
Customer	The customer raising the complaint	Required
Subject	Brief title of the complaint	Required
Description	Full details of the issue reported	Required
Status	Current resolution status	Required
Resolution Notes	How the complaint was resolved	

Customer Feedback

Customers can rate their experience via the Customer Portal on four dimensions. Feedback records appear in this module for review by the Quality Manager.

- Turnaround Time Rating (1–5 stars)
- Certificate Quality Rating (1–5 stars)
- Staff Behaviour Rating (1–5 stars)
- Overall Satisfaction Rating (1–5 stars)
- Comments (free text)

Lab Documents

Document control system for quality manual, SOPs, and records

ISO 17025 Compliance

Document Control

Maintain the lab's quality document library: Quality Manual, Standard Operating Procedures (SOPs), work instructions, forms, and records. The system manages version numbers and review dates.

Document Fields

Document No.	Unique document reference (e.g., QM-001, SOP-CAL-001)	Required
Title	Document title	Required
Category	Quality Manual / SOP / Work Instruction / Form / Record	Required
Version	Current version number (e.g., 1.0, 2.3)	Required
Revision Date	Date of last revision	Required
Next Review	Scheduled next review date — triggers notification when due	
File Attachment	Upload the document file (PDF, DOCX)	

Internal Audit

Planning and recording internal quality system audits

Internal Audit Module

Plan, execute, and record results of internal audits as required by ISO/IEC 17025. Link findings to NCR/CAPA records for follow-up.

Audit Fields

Audit Reference	Unique audit identifier	Required
Audit Date	Date the audit was conducted	Required
Auditor	Name of the internal auditor	Required
Area Audited	Department or scope (e.g., Calibration Lab, Document Control)	Required
Findings	Observations, non-conformances found, and opportunities for improvement	
Status	PLANNED → IN_PROGRESS → COMPLETED → CLOSED	

Tasks

Internal task management and assignment

| Tasks Module

Create and assign internal tasks to lab staff. Tasks are independent of calibration jobs — they cover any lab activity: maintenance, procurement follow-up, training, etc.

| Task Fields

Title	Task title	Required
Description	Detailed description of what needs to be done	
Assigned To	Select a user from the system	Required
Priority	LOW / MEDIUM / HIGH / URGENT	Required
Due Date	Deadline for completion	
Status	TODO → IN_PROGRESS → DONE → CANCELLED	Auto

CHAPTER 21

Engineers

Managing calibration engineers and their competency records

| Engineers Module

The Engineers module tracks all calibration and service engineers, their competencies, and links to their user account for job assignment.

| Engineer Fields

Employee Code	Unique staff identifier (e.g., ENG-001)	Required
Linked User Account	The system user account for this engineer	Required
Disciplines	Calibration disciplines this engineer is qualified for (Electrical, Mechanical, etc.)	
Qualification	Educational and professional qualifications	
Joining Date	Date of joining the lab	

Engineers must have a linked user account before they can be assigned to jobs. Create the user account first in Administration → Users.

CHAPTER 22

Users & Role Permissions

Creating staff accounts and configuring module-level access control

User Management

The LAB_ADMIN creates system users for all lab staff. Each user gets one of the six available roles.

Users — Administration + New User				
Name	Email	Role	Status	Actions
Dr. Anjali Kumar	anjali@lab.com	TECHNICAL_MANAGER	Active	  
Ravi Sharma	ravi@lab.com	CALIBRATION_ENGINEER	Active	  

Creating a User

Full Name	Staff member's full name	Required
Email	Login email address (must be unique)	Required
Role	One of: Technical Manager, Quality Manager, Calibration Engineer, Service Engineer, Data Entry Operator	Required
Password	Initial password (user should change on first login)	Required

Role Permissions

Go to **Administration** → **Role Permissions** to configure which modules each role can access. This is a grid of roles vs. modules with read/write toggles.

Role Permissions					
Module	Tech Mgr	Quality Mgr	Calib. Eng.	Svc. Eng.	Data Entry
Jobs	☒ Read/Write	☐ Read	☒ Read/Write	☒ Read/Write	☐ Read
Certificates	☒ Read/Write	☐ Read	☐ Read	☒ X	☒ X

Customers	✔ Read/Write	👁 Read	👁 Read	👁 Read	✔ Read/Write
NCR / CAPA	👁 Read	✔ Read/Write	✗	✗	✗
Complaints	👁 Read	✔ Read/Write	✗	✗	✗

Available modules in the permissions matrix:

- Jobs
- Certificates
- Instruments
- Reference Standards
- Calibration Masters
- Environmental
- Tasks
- Engineers
- Customers
- Quotations
- Invoices
- Purchase Orders
- Delivery Challans
- Leads / Pipeline
- CRM Activities
- NCR / CAPA
- Complaints & Feedback
- Lab Documents
- Internal Audit
- Reports
- Inventory
- Notifications

Inventory

Managing lab consumables and equipment stock

Inventory Module

Track quantities of consumables, chemicals, spare parts, and lab supplies. Set minimum stock levels to receive low-stock alerts.

Item Fields

Item Name	Name of the consumable or item	Required
Category	Consumable / Chemical / Spare Part / Equipment / Other	Required
Unit	Unit of measure (pcs, litres, kg, rolls)	Required
Current Qty	Current stock quantity	Required
Min. Stock Level	Quantity below which a low-stock alert is triggered	
Location	Storage location in the lab	

Items at or below minimum stock level are highlighted in orange and appear on the dashboard alerts panel.

Reports

Analytics, summaries, and exportable reports

Reports Module

Generate and export reports across all modules. Use date range filters, customer filters, and engineer filters to narrow down data.

Available Reports

Report	Description	Formats
Jobs Summary	Count of jobs by status, engineer, customer for a period	Table / CSV
Certificate Register	All certificates issued with validity dates	Table / CSV / PDF
Engineer Productivity	Number of jobs completed per engineer	Chart / CSV
Revenue Report	Invoice totals, GST collected, outstanding dues	Table / PDF
Customer-wise Revenue	Revenue breakdown per customer	Chart / CSV
NCR Summary	Non-conformances by category, status, and period	Table
Environmental Log	Temperature/humidity trends over time	Chart / CSV
Instrument Due List	Instruments with calibration due in a date range	Table / CSV

Notifications

System alerts and activity notifications

Notifications Panel

The bell icon in the top navigation bar shows unread notification count. Click to view all notifications.

Notification Types

Trigger	Who is Notified
New job received	Technical Manager, Admin
Job status changed	Assigned engineer, Technical Manager
Job pending review	Technical Manager
Certificate generated	Admin, customer (if email configured)
Reference standard expiring	Technical Manager, Admin
Instrument calibration due	Technical Manager, Admin
Low inventory alert	Admin
New complaint received	Quality Manager, Admin
NCR assigned to user	Assigned user
Internal audit due	Quality Manager, Admin

Admin Settings

Lab configuration, certificate templates, signatures, and data management

LAB_ADMIN Only

Accessing Settings

Navigate to **Administration** → **Settings**. This section is accessible only to LAB_ADMIN. It is organized into tabs:

Tab 1: Lab Details

Configure the lab's identity and contact information that appears on certificates, invoices, and quotations.

Lab Name	Official laboratory name	Required
Accreditation No.	NABL accreditation number	
Address	Full registered address	
Phone / Email	Lab contact information	
GSTIN / PAN	Tax registration numbers	
Logo	Click the logo box to upload. Appears on certificates and PDF documents. Accepted: PNG, JPG (max 2MB)	

Tab 2: Bank Details

Bank Name	Bank name (e.g., HDFC Bank)
Branch	Branch name
Account Number	Lab's bank account number
IFSC Code	Bank IFSC code

Tab 3: Certificate Signature Block

Configure the two signatures that appear at the bottom of every calibration certificate.

Technical Signatory (Calibration verifier)	Authorized Signatory (Approver / Director)
Technical Signatory Name e.g. Dr. A. Kumar	Authorized Signatory Name e.g. Mr. R. Sharma
Technical Signatory Designation e.g. Technical Manager	Authorized Signatory Designation e.g. Director
Signature Image Click to upload signature	Signature Image Click to upload signature

Signature Upload: Click the dashed box to open a file picker. Select a PNG or JPG of the handwritten signature. The image is stored as a data URL and printed on all certificates. Click *Remove* below the preview to clear it.

Tab 4: Certificate Templates

View and preview the available certificate templates. Click the **Preview** button on any template to see how the certificate looks with your lab's data.

Template	Used For	Includes NABL Logo
Template A — Standard NABL	Accredited calibrations	✓ Yes
Template B — Non-NABL	In-house calibrations	✗ No
Template C — Onsite	Onsite calibrations	Conditional
Template D — Multi-point	Multi-range calibrations	✓ Yes

Tab 5: Data Management

Export data from all modules as CSV for external reporting or backup.

- **Load Sample Data** — Populates the system with 5 records per module for demo/testing purposes
- **Export Customers** — Download all customer records as CSV
- **Export Instruments** — Download all instruments
- **Export Jobs** — Download all jobs with status
- **Export Certificates** — Certificate register as CSV
- **Export Invoices** — All invoices with amounts
- **Export Leads** — Sales pipeline export

Customer Portal

Self-service portal for customers to track jobs, download certificates, and raise complaints

What is the Customer Portal?

The Customer Portal is a separate login interface available at `/portal`. Customers use it to independently check their calibration job status, download certificates, submit complaints, and give feedback — without contacting the lab directly.

Portal Login

 Customer Portal — Login



Customer Portal Login

Track your calibration jobs

Lab Accreditation Code

Email Address

Password

Sign In

Setting Up Portal Access for a Customer

1. Go to **Customers** module
2. Find the customer and click the  (Lock) icon
3. Enter a portal password and click *Set Password*

4. Share the customer's registered email, the new password, and your lab's accreditation code with them
5. The customer logs in at /portal

Portal Tabs

Tab	What the Customer Sees / Can Do
 My Jobs	View all calibration jobs for their instruments: job number, instrument name, current status, date received
 Certificates	View all issued certificates. Click <i>Download</i> to get a PDF of the certificate with full datasheet
 My Instruments	View all instruments registered under their account with calibration due dates
 Raise Complaint	Submit a complaint with a subject and description. It instantly appears in the lab's <i>Quality</i> → <i>Complaints</i> module
 Feedback	Rate the lab on 4 dimensions (Turnaround Time, Certificate Quality, Staff Behaviour, Overall Satisfaction) and add comments

Portal Complaints: When a customer submits a complaint from the portal, it creates a record in the lab's Complaints module with the prefix PORTAL/ and is immediately visible to the Quality Manager and Lab Admin.

Certificate Download (Portal)

When a customer clicks Download on a certificate in the portal, the system fetches the full certificate data and generates a formatted HTML certificate that opens in a new browser tab, which can then be printed or saved as PDF.

APPENDIX

Quick Reference

Keyboard shortcuts, status codes, and glossary

Status Code Reference

Code	Module	Meaning
RECEIVED	Jobs	Instrument received at lab, job created
WAITING	Jobs	Waiting for engineer availability or parts
ASSIGNED	Jobs	Engineer has been assigned
IN_CALIBRATION	Jobs	Calibration work in progress
PENDING_REVIEW	Jobs	Engineer submitted; awaiting Technical Manager review
CORRECTION_REQUIRED	Jobs	Technical Manager flagged issues; engineer must re-do
APPROVED	Jobs	Technical Manager approved; ready for certificate
CERTIFICATE_GENERATED	Jobs	Certificate issued and available for download
DELIVERED	Jobs	Instrument returned to customer
CLOSED	Jobs	Job fully complete, no further action
NABL	Certificates	Issued under NABL accreditation scope
Non-NABL	Certificates	Issued outside NABL accreditation scope

Glossary

Term	Definition
NABL	National Accreditation Board for Testing and Calibration Laboratories — India's accreditation body
NCR	Non-Conformance Report — document describing a deviation from defined quality requirements
CAPA	Corrective and Preventive Action — actions taken to eliminate root causes and prevent recurrence
DC	Delivery Challan — goods movement document
Measurand	The physical quantity being measured during calibration
Traceability	Chain linking calibration results back to national/international measurement standards

MUT	Measurement Under Test — the instrument being calibrated
CMC	Calibration and Measurement Capability — the lab's best achievable uncertainty
SOP	Standard Operating Procedure
Portal	The customer-facing self-service web interface at /portal

Support

For technical support or feature requests, contact your system administrator or the CLMS support team.



CortexCLMS

Complete User Manual · Version 1.0 · 2026

This document covers all modules and functionality available in CLMS v1.0

 +91 8329925318

 cortexaitechnologies@zohomail.in

 Pune, Maharashtra, India 411034